

**General guidelines for laboratory
detection of carbapenemase-producing
Enterobacterales (CPE) in Belgium**

NRC Antibiotic-resistant Gram-negative bacilli

Carbapenemase types in *Enterobacterales*, *Pseudomonas* spp. and *Acinetobacter* spp. vs

Amblar molecular class	Carbapenemase family	<i>Enterobacterales</i> (CPE)	<i>Pseudomonas</i> spp.	<i>Acinetobacter</i> spp.
A (serine)	KPC GES IMI, NMC, SME, SFC...	+++ (+) +	- + -	- (+) -
B (zinc MBL)	VIM, NDM IMP AIM, DIM, GIM, SIM, SPM...	+++ +++ + -	+++ + + (+)	+ ++ + (+)
D (serine)	OXA-48-like OXA-427 OXA-198 OXA-23-like OXA-24-like OXA-58-like	+++ + (+) (+) (+) (+)	- - + - - -	- - - +++ ++ ++

Selection criteria for suspicion of CPE

- Species identification confirmed belonging to *Enterobacterales*
- Decreased susceptibility to at least one carbapenem using screening breakpoint (not clinical S breakpoint):
 - MIC method: **ertapenem OR meropenem > 0.125 mg/L** (Sensitivity of 95%)
 - NB: automated systems may have insufficient sensitivity if lowest concentration does not reach 0.125 mg/l
 - Disk diffusion method:
 - **Ertapenem (10 µg) < 25 mm** (Sensitivity of 98-99%) OR
 - **Meropenem (10 µg) < 28 mm** (Sensitivity of 95%)

If carba non-S confirmed → Perform first-line screening tests

First-line screening tests to be performed by local laboratories

At least one test positive of the following :

- **High-level resistance to ceftazidime/avibactam CAZ/AVI :**
 - Disk diffusion: CAZ/AVI (10/4 µg) <10 mm OR
 - MIC method: CAZ/AVI >16 mg/L

OR

- **High-level resistance to temocillin TMO :**
 - Disk diffusion: TMO (30 µg) <12 mm OR
 - MIC method: TMO >64 mg/L

OR

- **Carbapenemase hydrolysis-based test (HBT) : one of the following**
 - **Colorimetric tests:**
 - Rapid CARB Screen Kit, Rapid CARB Blue Kit (ROSCO)
 - RAPIDEC (BioMérieux)
 - Beta-CARBA test (BioRad)
 - **Mass spectrometry tests:**
 - MALDITOF MBT Star Carba (Bruker)
 - **Carbapenem inactivation methods (CIM)**
 - Modified CIM : mCIM, zCIM...

Second-line confirmatory tests (optional)

1) Immunochromatographic tests (ICT):

- Monoplex ICT: OXA-48 → 60% of all CPE in Belgium
- Multiplex ICT: OXA-48, KPC, NDM, VIM, IMP

2) Inhibitor-based combination (IBT) disk tests or gradient MIC tests:

- Class B carbapenemase (MBL: VIM/NDM/IMP): Mero +/- DPA (dipicolinic acid) or Imi +/- EDTA
- Class A carbapenemase (KPC): Mero +/- BA (boronic acid) and mero +/- cloxacillin

3) Molecular tests (PCR assays, DNA microarray, isothermal amplification, sequencing...): identification of carbapenemase family and/or enzyme variant

General criteria to rejection of request or delay of analysis for confirmation of CPE by the NRC

- Sample received without reference ID, minimal sample **information** (sample nature / date), isolate species identification, copy of antibiogram results,...
- **Damaged** samples (crushed plates, broken glass tubes,....)
- **Contaminated** culture
- Isolate received at the NRC identified as a **different species** from the one notified on the request form
- Isolate belonging to genera/species **other than *Enterobacterales***
- Isolate with **no decreased susceptibility** to any carbapenem (ertapenem, meropenem) according to screening criteria
- **No first-line nor second-line tests** being performed by local laboratories
- Repeated isolates of **same species from the same patient** with similar resistance patterns

Algorithm for detection of CPE on isolates

